

## Product data sheet

29 April 2020

# 1905 Isolation Gown

Isolation Gown

### Product details

Size : L  
Colour : Dark Blue  
Sterile : Non-sterile

### Images



### Delivered items

1905-02

Country of origin: Cambodia

Shelf life: 5 years

Sterilization method: Non-sterile

Packing information: First packaging layer is a plastic liner/bag. Second layer is a corrugated board transport box. Is suitable for Tray: No

Packing level	Quantity	GS1 code	WxLxH (mm)	Vol (dm3)	Weight gross/net (kg)
Transport box	50	7323190050055	305x330x330	33.2	3.4 / 2.9
Pallet	1500	7323190050048	800x1200x1675		

### Material

Animal tissues : No  
Natural rubber latex : No  
Medicinal substances : No

### Product Composition Wearing Apparel

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
©2020 Mölnlycke Health Care AB. All rights reserved.



## Product data sheet

29 April 2020

Product Component	Composition
Main material	Polypropylene spunbond nonwoven 28 g/m <sup>2</sup>
Product Component	Composition
Neck band	Polypropylene spunbond nonwoven 28 g/m <sup>2</sup>
Elastic	Polyurethane, Polyamide (Nylon)
Tie bands	Polypropylene spunbond nonwoven 28 g/m <sup>2</sup>
Sewing thread	Polyurethane, Polyamide (Nylon)

### Technical

#### Dimension

Dimension text	Dimension value
Length, highest point at shoulder	114 cm
Center front	103 cm
Chest width	65 cm
Top sleeve	56 cm
Armhole width	28 cm

#### Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

#### Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

#### Classifications

Regulation type(s)	MDD Class I
MDD Classification Rule :	1
Intended use MDD :	The intended use is to reduce contamination in a non-sterile environment within health care.

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.  
The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.  
©2020 Mölnlycke Health Care AB. All rights reserved.



## Product data sheet

29 April 2020

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 10993-5, EN ISO 15223-1, ISO 15223-2, EN ISO 10993-10, ISO 14001

Removable label

No

GMDN Code (Global Medical Device Nomenclature)

35492 Isolation gown, single-use

UNSPSC

42131612 Medical staff isolation or cover gowns

Commodity Code

6210109200 Wearing apparel mainly made of nonwovens

Find out more at [www.molnlycke.com](http://www.molnlycke.com)

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.

The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2020 Mölnlycke Health Care AB. All rights reserved.

